


AS 50-542 Retractor Kit II

Includes

- AS 50-552 Horizontal Bar
- AS 50-554 1" Post Coupling
- AS 50-557 Medium Oval Ring
- AS 50-561 2" x 6" Malleable Retractors
- AS 50-562 11½ x 6" Malleable Retractors 38 x 152 mm
- AS 50-566 Harrington Retractor
- AS 50-567 3" x 6" Malleable Retractors 76 x 152 mm
- AS 50-568 2" x 3" Kelly Retractors (2ea) 51 x 76 mm
- AS 50-569 2" x 4" Kelly Retractors 51 x 101 mm
- AS 50-570 2" x 5" Kelly Retractors 51 x 127 mm
- AS 50-571 2" x 6" Kelly Retractors 51 x 152 mm
- AS 50-576 Balfour Retractor
- AS 50-577 Gelpi Retractor
- AS 50-579 Ratchet Mechanismus (2ea)
- AS 50-580 Tilt Ratchet Mechanismus (6ea)
- AS 50-582 Horizontal Flex Bar
- AS 50-583 Segmented Ring medium I.D.
10½" x 19½" (266 x 495 mm)



Optional

- AS 50-547 Rectal Kit
- AS 50-553 Extension Part (Cross pole)
- AS 50-556 Oval Rings large width(I.D) 12"
(267 mm) length (I.D) 18½" (470 mm)
- AS 50-558 Oval Rings small width (I.D) 8" (230 mm)
length (I.D) 12¼" (311 mm)
- AS 50-559 Round Rings medium 10½" I.D (266 mm)
- AS 50-560 Round Rings large 12½" I.D (317 mm)
- AS 50-563 Kelly Retractor Blade 51 x 51 mm
- AS 50-565 Kelly Blade Modified width 2" (51 mm)
depth 1½" (38 mm)
- AS 50-572 Kelly Blade Standart width 2" (51 mm)
depth 2½" (64 mm)
- AS 50-574 Kelly Blade 1" x 1½" (25 x 38 mm)
- AS 50-575 Kelly Blade 1" x 1½" (25 x 63 mm)
- AS 50-581 Table Post
- AS 50-584 Round Ring small 8½" I.D. (216 mm)
- AS 50-585 Rectal Blades small (2ea) 3" x 7" (76 x 78 mm)
- AS 50-586 Rectal Blades large (2ea) 4" x 6" (101 x 152 mm)
- AS 50-587 Segmented Ring small I.D. 8½" x 14" (216 x 356 mm)
- AS 50-588 Segmented Ring large 12½" x 19½" (317 x 495 mm)
- AS 50-590 Vaginal Retractor Ring I.D. 8" x 10½" (203 x 206 mm)
- AS 50-591 Lateral Retractor Blade (2ea) 1" x 3" (25 x 76 mm)
- AS 50-592 Deaver Blade small 1" x 3½" (25 x 89 mm)
- AS 50-593 Deaver Blade large 1" x 5" (25 x 127 mm)
- AS 50-594 Posterlor Blade 1" x 4½" (114 mm)
- AS 50-635 Segment Ring round Ø216 mm
- AS 50-639 Scoville Hook Retractor small (2each)
- AS 50-640 Scoville Hook Retractor large (2each)
- AS 50-641 Meyerding Retractor small (2each)
- AS 50-642 Meyerding Retractor medium (2each)
- AS 50-643 Meyerding Retractor large (2each)
- AS 50-644 Deep Rake Retractor small (2each)
- AS 50-645 Skin Rake Retractor medium (2each)
- AS 50-646 Skin Rake Retractor large (2each)

